

New technology for skin cancer prevention





Mole Mapping versus Total Body Mapping



Mole MappingSingle moles or body areas



Total Body Mapping

➤ Complete skin surface





- > Documentation of the whole skin surface
- Reproducible photos from head to toe
- > Standardized and systematic imaging process
- > Comparable results





- > Documentation of the whole skin surface
- > Reproducible photos from head to toe
- > Standardized and systematic imaging process
- Comparable results





- > Documentation of the whole skin surface
- > Reproducible photos from head to toe
- > Standardized and systematic imaging process
- > Comparable results





- > Documentation of the whole skin surface
- > Reproducible photos from head to toe
- > Standardized and systematic imaging process
- > Comparable results




New: FotoFinder ATBM® - Automated Total Body Mapping





The 4 pillars of ATBM® skin cancer prevention





FotoFinder Automated Total Body Mapping



The device...



The new FotoFinder bodystudio ATBM®



- ➤ German innovation in skin cancer diagnosis
- ➤ 1st system for Automated Total Body Mapping in the world
- Seamless integration of dermoscopy



A revolution in skin cancer diagnosis!



- ➤ Total Body Mapping in less than 3 minutes
- > Fast detection of skin changes
- Intuitive software-guided procedure that can be delegated



ATBM[®] camera system



- > Fully automated camera positioning
- ➤ Software-controlled digital SLR camera
- ➤ Unique BodyFlash system
- SmartMatch technology for whole body views with up to 72 megapixels



Full HD camera for dermoscopy



- > Full HD live video camera
- Variable optical magnification
- > Autofocus for all zoom levels
- > Immersion or polarization
- ➤ Unachieved image quality



Reproducible photos at baseline and follow-ups



- ➤ Laser line: reproducible distance
- > Floor mat: easy foot positions
- ➤ Ghost: reproducible patient pose



FotoFinder quality hardware



- ➤ Silent Medical Server
- Multiple backup system
- ➤ 23" quality screen
- ➤ Mobile tower with stable castors




FotoFinder Automated Total Body Mapping



The procedure...



20 standard positions plus additional segments






































The procedure takes only 3 minutes!






Worldwide unique: ATBM® capturing procedure



- > The procedure can be done by your staff
- ➤ It takes 20 seconds per body view (4 photos)
- > After 3 minutes....
 - o 20 photos are taken
 - Complete Bodyscan analysis is done for all 20 pictures
 - All information is available in the network on several workstations




New Bodyscan ATBM® shows changes






Details visible already in the total body photo






Seamless integration of dermoscopy





Mole history at a glance






Brilliant microscopic images in Full HD live quality






Brilliant microscopic images in Full HD live quality






Brilliant microscopic images in Full HD live quality





Moleanalyzer score: high sensitivity and specificity





ATBM® – Advantages at a glance



- > Total body views with 72 megapixels
- > Fully automated camera movement
- ➤ Laser positioning system
- > Extremely fast indication of skin changes
- ➤ Unrivalled Full HD quality dermoscopy
- ➤ Mole analysis: high sensitivity <u>and</u> specificity
- Comprehensive patient reports
- > Easy procedure can be done by the nurse
- Future-proof FotoFinder technology
- Imaging quality made in Germany



FotoFinder. Setting the Pace in Skin Imaging Worldwide!



- > Pioneer since more than 20 years
- World leading technology for video dermoscopy and total body mapping
- > Inventor of eDermoscopy and handyscope
- Several thousand customers in more than 70 countries.
- > Family-owned business
- Certified quality management according to ISO 9001:2008